

Product Datasheet

HA peptide (orb216634)

Description	HA peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Tyr-Pro-Tyr-Asp-Val-Pro-Asp-Tyr-Ala-OH
Formula	$C_{53}H_{67}N_9O_{17}$
Purity	> 95%
Source	Synthetic
MW	1102.18
CAS Number	92000-76-5
Expiration Date	12 months from date of receipt.

