

Product Datasheet

L-Selectin peptide (orb216676)

Description	L-Selectin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Cys-Gln-Lys-Leu-Asp-Lys-Ser-Phe-Ser-Met-Ile-Lys-OH
Formula	$C_{62}H_{106}N_{16}O_{18}S_2$
Purity	> 95%
Source	Synthetic
MW	1427.73
CAS Number	126880-86-2
Expiration Date	12 months from date of receipt.