

Product Datasheet

Motilin, Canine peptide (orb216707)

Description	Motilin, Canine peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	A gastric motor activity-stimulating polypeptide.
Protein Sequence	H-Phe-Val-Pro-Ile-Phe-Thr-His-Ser-Glu-Leu-Gln-Lys-Ile-Arg-Glu-Lys-Glu-Arg-Asn-Lys-Gly-Gln-OH
Formula	$C_{120}H_{194}N_{36}O_{34}$
Purity	> 95%
Source	Synthetic
MW	2685.1
CAS Number	85490-53-5
Expiration Date	12 months from date of receipt.

