

Product Datasheet

OVA (323-339) peptide (orb216747)

Description	OVA (323-339) peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	ISQAVHAAHAEINEAGR is used for studies of T-cell activation and allergic diseases.
Protein Sequence	Ile-Ser-Gln-Ala-Val-His-Ala-Ala-His-Ala-Glu-Ile-Asn-Glu-Ala-Gly-Arg
Formula	$C_{74}H_{120}N_{26}O_{25}$
Purity	> 95%
Source	Synthetic
MW	1773.9
CAS Number	92915-79-2
Expiration Date	12 months from date of receipt.

