

Product Datasheet

P55-TNFR peptide (orb216750)

Description	P55-TNFR peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Leu-Pro-Gln-Ile-Glu-Asn-Val-Lys-Gly-Thr-Glu-Asp-OH
Formula	$C_{57}H_{95}N_{15}O_{22}$
Purity	> 95%
Source	Synthetic
MW	1342.48
Expiration Date	12 months from date of receipt.