

Product Datasheet

Prolactin peptide (orb216771)

Description	Prolactin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Prolactin-releasing peptide. Stimulates luteinizing hormone and follicle stimulating hormone.
Protein Sequence	Ser-Arg-Thr-His-Arg-His-Ser-Met-Glu-Ile-Arg-Thr-Pro-Asp-Ile-Asn-Pro-Ala-Trp-Tyr-Ala-Ser-Arg-Gly-Ile-Arg-Pro-Val-Gly-Arg-Phe-NH ₂
Formula	C ₁₆₀ H ₂₅₂ N ₅₆ O ₄₂ S
Purity	> 95%
Source	Synthetic
MW	3664.2
CAS Number	215510-22-8
Expiration Date	12 months from date of receipt.