

Product Datasheet

Ranatensin R peptide (orb216778)

Description	Ranatensin R peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Asn-Thr-Ala-Leu-Arg-Arg-Tyr-Asn-Gln-Trp-Ala-Thr-Gly-His-Phe-Met-NH ₂
Formula	C ₉₀ H ₁₃₄ N ₃₀ O ₂₄ S
Purity	> 95%
Source	Synthetic
MW	2052.3
CAS Number	70572-93-9
Expiration Date	12 months from date of receipt.

