

Product Datasheet

Renin Inhibitor peptide (orb216779)

Description	Renin Inhibitor peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-His-Pro-Phe-His-Leu-D-Leu-Val-Tyr-NH ₂
Formula	C ₅₂ H ₇₃ N ₁₃ O ₉
Purity	> 95%
Source	Synthetic
MW	1024.24
Expiration Date	12 months from date of receipt.

