

Product Datasheet

R-G-D-S peptide (orb216786)

Description	R-G-D-S peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	RGDS, a tetrapeptide which supports fibroblast attachment and inhibits fibronectin binding to platelets. RGDS is a target sequence for syphilis spirochete cytodherence.
Protein Sequence	Arg-Gly-Asp-Ser
Formula	$C_{15}H_{27}N_7O_8$
Purity	> 95%
Source	Synthetic
MW	433.4
CAS Number	91037-65-9
Expiration Date	12 months from date of receipt.

