

Product Datasheet

R-S-R peptide (orb216792)

Description	R-S-R peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Tholander and Haeggström developed an assay for the rapid determination of the tripeptidase activity of LTA4 hydrolase using RSR.
Protein Sequence	H-Arg-Ser-Arg-OH
Formula	$C_{15}H_{31}N_9O_5$
Purity	> 95%
Source	Synthetic
MW	417.5
CAS Number	115035-42-2
Expiration Date	12 months from date of receipt.

