

Product Datasheet

TAT 2-4 peptide (orb216829)

Description	TAT 2-4 peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Tyr-Gly-Arg-Lys-Lys-Arg-Arg-Gln-Arg-Arg-Arg-Gly-Tyr-Gly-Arg-Lys-Lys-Arg-Arg-Gln-Arg-Arg-Arg-Gly-OH
Formula	$C_{132}H_{240}N_{66}O_{29}$
Purity	> 95%
Source	Synthetic
MW	3215.81
Expiration Date	12 months from date of receipt.

