

Product Datasheet

Triptorelin peptide (orb216837)

Description	Triptorelin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	LHRH agonist used in the treatment of prostate cancer. Triptorelin decreases the secretion of LH and FSH.
Protein Sequence	pGlu-His-Trp-Ser-Tyr-D-Trp-Leu-Arg-Pro-Gly-NH ₂
Formula	C ₆₄ H ₈₂ N ₁₈ O ₁₃
Purity	> 95%
Source	Synthetic
MW	1311.46
CAS Number	57773-63-4
Expiration Date	12 months from date of receipt.

