

Product Datasheet

Universal TT epitope P2 (830-844) peptide (orb216840)

Description	Universal TT epitope P2 (830-844) peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gln-Tyr-Ile-Lys-Ala-Asn-Ser-Lys-Phe-Ile-Gly-Ile-Thr-Glu-Leu
Formula	$C_{80}H_{129}N_{19}O_{23}$
Purity	> 95%
Source	Synthetic
MW	11725.03
Expiration Date	12 months from date of receipt.