

Product Datasheet

Uperolein peptide (orb216841)

Description	Uperolein peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	pGlu-Pro-Asp-Pro-Asn-Ala-Phe-Tyr-Gly-Leu-Met-NH ₂
Formula	C ₅₇ H ₇₉ N ₁₃ O ₁₆ S
Purity	> 95%
Source	Synthetic
MW	1234.42
Expiration Date	12 months from date of receipt.

