

Product Datasheet

Urocortin III, Mouse peptide (orb216845)

Description	Urocortin III, Mouse peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	High affinity, selective CRF2 receptor agonist. Stimulates glucagon and insulin release from isolated rat islets.
Protein Sequence	H-Phe-Thr-Leu-Ser-Leu-Asp-Val-Pro-Thr-Asn-Ile-Met-Asn-Ile-Leu-Phe-Asn-Ile-Asp-Lys-Ala-Lys-Asn-Leu-Arg-Ala-Lys-Ala-Ala-Ala-Asn-Ala-Gln-Leu-Met-Ala-Gln-Ile-NH ₂
Formula	C ₁₈₆ H ₃₁₂ N ₅₂ O ₅₂ S ₂
Purity	> 95%
Source	Synthetic
MW	4173.01
CAS Number	357952-10-4
Expiration Date	12 months from date of receipt.