

Product Datasheet

Xenin peptide (orb216859)

Description	Xenin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Xenin 25, which has 6 C-terminal amino acids in common with amphibian xenopsin, has been detected in human gastric mucosa. It stimulates the exocrine pancreatic secretion.
Protein Sequence	H-Met-Leu-Thr-Lys-Phe-Glu-Thr-Lys-Ser-Ala-Arg-Val-Lys-Gly-Leu-Ser-Phe-His-Pro-Lys-Arg-Pro-Trp-Ile-Leu-OH
Formula	$C_{139}H_{224}N_{38}O_{32}S$
Purity	> 95%
Source	Synthetic
MW	2971.63
CAS Number	144092-28-4
Expiration Date	12 months from date of receipt.

