

Product Datasheet

A1BG peptide (orb216864)

Description	The peptide is used to block anti-Alpha 1B glycoprotein antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	A1BG
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213508.
Purity	> 85%
Source	Synthetic
Uniprot ID	P04217 , Q9EPH1 , Q19LI2
Entrez	140656 , 117586 , 1
Expiration Date	12 months from date of receipt.