

Product Datasheet

ABL1 peptide (orb216868)

Description	The peptide is used to block anti-ABL1/2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ABL1; ABL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213512.
Purity	> 85%
Source	Synthetic
Uniprot ID	P42684 , Q4JIM5 , P00520 , P00519
Entrez	27 , 11350 , 25
Expiration Date	12 months from date of receipt.