

Product Datasheet

ACAT1 peptide (orb216869)

Description	The peptide is used to block anti-ACAT1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	ACAT1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213513.
Purity	> 85%
Source	Synthetic
Uniprot ID	P24752, P17764, Q8QZT1
Entrez	110446, 25014, 38
Expiration Date	12 months from date of receipt.