

Product Datasheet

ACHE peptide (orb216870)

Description	The peptide is used to block anti-Acetylcholinesterase antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	ACHE
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213514.
Purity	> 85%
Source	Synthetic
Uniprot ID	P22303, P21836, P37136
Entrez	43, 83817, 11423
Expiration Date	12 months from date of receipt.