

Product Datasheet

beta Actin peptide (orb216872)

Description	beta Actin peptide
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	ACTB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213516.
Purity	> 85%
Source	Synthetic
Uniprot ID	P60710, P60711, P60709
Entrez	11461, 81822, 60
Expiration Date	12 months from date of receipt.