

Product Datasheet

ACTN3 peptide (orb216876)

Description	The peptide is used to block anti-ACTN3 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	ACTN3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213520.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q08043, O88990
Entrez	11474, 89
Expiration Date	12 months from date of receipt.