

Product Datasheet

ACVR1B peptide (orb216877)

Description	The peptide is used to block anti-ACVR1B antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ACVR1B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213521.
Purity	> 85%
Source	Synthetic
Uniprot ID	P36896, Q61271, P80202
Entrez	29381, 91, 11479
Expiration Date	12 months from date of receipt.