

Product Datasheet

ADCY5 peptide (orb216882)

Description	The peptide is used to block anti-Adenylate Cyclase 5/6 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	ADCY5; ADCY6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213526.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q01341 , O43306 , O95622 , Q03343 , P84309 , Q04400
Entrez	11512 , 224129 , 111 , 64532 , 112
Expiration Date	12 months from date of receipt.