

Product Datasheet

ADD3 peptide (orb216884)

Description	The peptide is used to block anti-Gamma-adducin antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ADD3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213528.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UEY8, Q62847, Q9QYB5
Entrez	25230, 120, 27360
Expiration Date	12 months from date of receipt.