

Product Datasheet

ADH7 peptide (orb216885)

Description	The peptide is used to block anti-ADH7 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	ADH7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213529.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q64437, P40394, P41682
Entrez	171178, 131, 11529
Expiration Date	12 months from date of receipt.