

Product Datasheet

ADORA2A peptide (orb216886)

Description	ADORA2A peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	ADORA2A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213530.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q60613, P29274, P30543
Entrez	135, 25369, 11540
Expiration Date	12 months from date of receipt.