

Product Datasheet

PARP1 peptide (orb216887)

Description	The peptide is used to block anti-PARP1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PARP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213531.
Purity	> 85%
Source	Synthetic
Uniprot ID	P11103, P27008, P09874
Entrez	142, 25591
Expiration Date	12 months from date of receipt.