

Product Datasheet

ADRA1B peptide (orb216889)

Description	The peptide is used to block anti-Alpha-1B Adrenergic Receptor antibody reactivity.
Reactivity	Human, Rat
Tested Applications	Blocking
Target	ADRA1B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213533.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15823, P35368
Entrez	147
Expiration Date	12 months from date of receipt.