

Product Datasheet

ADRBK1 peptide (orb216892)

Description	The peptide is used to block anti-GRK2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	ADRBK1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213536.
Purity	> 85%
Source	Synthetic
Uniprot ID	P26817, P25098, Q99MK8
Entrez	25238, 156
Expiration Date	12 months from date of receipt.