

Product Datasheet

ACAN peptide (orb216893)

Description	The peptide is used to block anti-Aggrecan antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	ACAN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213537.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61282, P16112, P07897
Entrez	58968, 11595
Expiration Date	12 months from date of receipt.