

Product Datasheet

AKT peptide (orb216894)

Description	AKT peptide
Reactivity	Human, Mouse, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	AKT1; AKT2; AKT3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213538.
Purity	> 85%
Source	Synthetic
Uniprot ID	P31750 , Q63484 , P47197 , P31751 , P47196 , Q9Y243 , Q60823 , Q9WUA6 , P31749
Entrez	23797 , 11651 , 207 , 29414 , 25233 , 10000 , 24185 , 208 , 11652
Expiration Date	12 months from date of receipt.