

Product Datasheet

AKT peptide (orb216897)

Description	AKT peptide
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	AKT1; AKT2; AKT3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213541.
Purity	> 85%
Source	Synthetic
Uniprot ID	P31750, Q63484, P47197, P31751, P47196, Q9Y243, Q60823, Q9WUA6, P31749
Entrez	23797, 11651, 207, 29414, 25233, 10000, 24185, 208, 11652
Expiration Date	12 months from date of receipt.