

Product Datasheet

AKT peptide (orb216901)

Description	AKT peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	AKT1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213545.
Purity	> 85%
Source	Synthetic
Uniprot ID	P47196, P31750, P31749
Entrez	207, 24185, 11651
Expiration Date	12 months from date of receipt.