

Product Datasheet

ALDOA peptide (orb216907)

Description	The peptide is used to block anti-ALDOA antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	ALDOA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213551.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05064, P05065, P04075
Entrez	11674, 226, 24189
Expiration Date	12 months from date of receipt.