

Product Datasheet

CD13 peptide (orb216910)

Description	CD13 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ANPEP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213554.
Purity	> 85%
Source	Synthetic
Uniprot ID	P97449, P15144
Entrez	16790, 290
Expiration Date	12 months from date of receipt.