

Product Datasheet

ANXA3 peptide (orb216914)

Description	The peptide is used to block anti-Annexin A3 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	ANXA3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213558.
Purity	> 85%
Source	Synthetic
Uniprot ID	P14669, O35639, P12429
Entrez	11745, 306, 25291
Expiration Date	12 months from date of receipt.