

Product Datasheet

AOX1 peptide (orb216916)

Description	The peptide is used to block anti-Aldehyde Oxidase antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	AOX1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213560.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q06278, Q9Z0U5, O54754
Entrez	11761, 316, 54349
Expiration Date	12 months from date of receipt.