

Product Datasheet

APAF1 peptide (orb216917)

Description	The peptide is used to block anti-APAF-1 antibody reactivity.
Reactivity	Human, Rat
Tested Applications	Blocking
Target	APAF1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213561.
Purity	> 85%
Source	Synthetic
Uniprot ID	O14727
Entrez	317
Expiration Date	12 months from date of receipt.