

Product Datasheet

APC peptide (orb216920)

Description	The peptide is used to block anti-APC antibody reactivity.
Reactivity	Human, Porcine, Rat
Tested Applications	Blocking
Target	APC
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213564.
Purity	> 85%
Source	Synthetic
Uniprot ID	P25054, P70478
Entrez	324, 24205
Expiration Date	12 months from date of receipt.