

Product Datasheet

APEX1 peptide (orb216921)

Description	The peptide is used to block anti-APE1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	APEX1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213565.
Purity	> 85%
Source	Synthetic
Uniprot ID	P28352, P43138, P27695
Entrez	79116, 11792, 328
Expiration Date	12 months from date of receipt.