

## Product Datasheet

### Apolipoprotein E peptide (orb216927)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Apolipoprotein E peptide                                 |
| <b>Reactivity</b>          | Human, Mouse, Primate, Rat                               |
| <b>Tested Applications</b> | Blocking                                                 |
| <b>Target</b>              | APOE                                                     |
| <b>Form/Appearance</b>     | Lyophilized powder                                       |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.             |
| <b>Note</b>                | For research use only                                    |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb213571. |
| <b>Purity</b>              | > 85%                                                    |
| <b>Source</b>              | Synthetic                                                |
| <b>Uniprot ID</b>          | <b>P08226, P02650, P02649</b>                            |
| <b>Entrez</b>              | <b>11816, 348, 25728</b>                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                          |