

Product Datasheet

Apolipoprotein E peptide (orb216928)

Description	Apolipoprotein E peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	APOE
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213572.
Purity	> 85%
Source	Synthetic
Uniprot ID	P08226, P02650, P02649
Entrez	11816, 348, 25728
Expiration Date	12 months from date of receipt.