

Product Datasheet

Aquaporin 4 Peptides (orb216932)

Description	Aquaporin 4 Peptides
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	AQP4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213576.
Purity	> 85%
Source	Synthetic
Uniprot ID	P55087, P55088, P47863
Entrez	25293, 11829, 361
Expiration Date	12 months from date of receipt.