

Product Datasheet

AR peptide (orb216934)

Description	The peptide is used to block anti-Androgen Receptor
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	AR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213578.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10275, P15207, P19091
Entrez	11835, 367, 24208
Expiration Date	12 months from date of receipt.