

Product Datasheet

ARAF peptide (orb216935)

Description	The peptide is used to block anti-A-RAF
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	ARAF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213579.
Purity	> 85%
Source	Synthetic
Uniprot ID	P04627, P10398, P14056
Entrez	369, 64363, 11836
Expiration Date	12 months from date of receipt.