

Product Datasheet

ARF6 peptide (orb216936)

Description	The peptide is used to block anti-ARF6 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	ARF6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213580.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62331, P62332, P62330
Entrez	11845, 382, 79121
Expiration Date	12 months from date of receipt.