

Product Datasheet

ARG1 peptide (orb216937)

Description	The peptide is used to block anti-Arginase 1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ARG1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213581.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05089, P07824, Q61176
Entrez	383, 29221, 11846
Expiration Date	12 months from date of receipt.