

Product Datasheet

ARG2 peptide (orb216938)

Description	The peptide is used to block anti-Arginase 2 antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	ARG2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213582.
Purity	> 85%
Source	Synthetic
Uniprot ID	P78540
Entrez	384
Expiration Date	12 months from date of receipt.