

Product Datasheet

ARSA peptide (orb216942)

Description	The peptide is used to block anti-Arylsulfatase A antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ARSA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213586.
Purity	> 85%
Source	Synthetic
Uniprot ID	P50428, P15289
Entrez	11883, 410
Expiration Date	12 months from date of receipt.